**Thera**vital

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**Instructions for use 35 ... 68**

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











Blumenweg 8  
D-88454 Hochdorf  
Tel.: +49 (0)7355-9314-0  
Fax: +49 (0)7355-9314-15  
e-mail: [info@thera-trainer.de](mailto:info@thera-trainer.de)  
<http://www.thera-trainer.de>






80 000 066-006 (Software 01.04.00-08/10/06)

**CE**  
0297

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Congratulations!

You have made an excellent choice with the purchase of the THERA-vital. This therapeutic exerciser is an innovative, high-quality product made in Germany. The THERA-vital is the first therapeutic exerciser to feature a colour display screen. It also functions according to the very latest therapeutic advances.

## Contents

Exerciser Options .....	36
Intended Use .....	36
Exerciser Components and Controls .....	37
Interfaces .....	37
Safety .....	38
Symbols Used .....	39
Notes on Exercise Therapy .....	40
General observations.....	40
Uses.....	40
Correct posture.....	40
Exercise programme planning .....	40
Active and active-assisted exercising ...	41
Muscle hypertonia/spasms .....	41
Transporting and Commissioning .....	42
Before using for the first time .....	42
Transporting.....	42
Connecting to the power supply .....	42
Switching the exerciser on/off .....	42
Operation.....	42
Foot rests.....	42
Radius adjustment 2-step .....	42

These instructions will help you to familiarise yourself with the THERA-vital. They will guide you carefully through its functions and operation and offer numerous tips and hints on getting the best use out of your new therapeutic exerciser.

Before using the exerciser for the first time, please take careful note of the safety instructions on page 38.

Radius adjustment stepless (optional feature).....	43
Tip-up protection (optional feature).....	43
Tip-up protection, self operated (optional feature).....	43
Calf rests (optional feature) .....	43
Foot fixings (optional feature) .....	43
Adjusting the handle bar .....	44
Adjusting the apparatus foot .....	44
Upper torso trainer (optional feature)....	44
Therapy grips for upper torso trainer ....	45
Safety handle adapters.....	45
Adjusting the position of the control pad..	45
Fore arm rests (optional feature) .....	46
Wristband (optional feature) .....	46
Heart rate sensor.....	46
Control pad with colour display.....	47
THERA-assist exercising software .....	49
Submenu (for experienced users) .....	50
Navigating the submenu .....	50
Menu structure .....	51
Description of exerciser functions on submenu.....	52

If you have any questions or observations, the staff at medica Medizintechnik GmbH will be only too pleased to help.

Enjoy exercising and stay active with the THERA-vital.

Exercising time .....	52
Forward/backward automatic.....	52
Biofeedback leg exerciser.....	52
Biofeedback arm exerciser .....	54
Motor power .....	55
Spasm detection .....	55
Spasm program .....	56
Automatic speed adaption .....	56
Pulse limit .....	56
Muscle support .....	56
Microphone sensitivity .....	56
Contrast .....	56
Keypad lock .....	57
Analysis of Session Data.....	57
Activating Basic Exercising Settings.....	58
Technical Data .....	58
Cleaning and Disinfection.....	59
Replacing a fuse .....	59
Guarantee Terms .....	59
Recycling .....	59
EC Declaration of Conformity .....	59
Tips and Tricks .....	60
Notes on electromagnetic compatibility .....	62
Index .....	66

## Exerciser Options

**Basic exerciser:** The basic version of the THERA-vital can be used for active-assisted (driven by the exerciser) and active (where the motive power is provided by the patient) therapeutic exercise of the legs. The patient's legs are held in the foot rests by Velcro straps. There is also an adjustable safety handle for balance and stability. All important parameters can be set on the control pad. The training data is displayed on the large colour display. The cranks can be set to two different lengths with the aid of suitable tools.

**Upper torso trainer option:** Instead of the safety handle, the THERA-vital can be fitted with an upper torso trainer. This has a separate motor, in other words both active-assisted and active exercise are possible. It is controlled from the same control pad as on the basic exerciser. The patient holds onto the exercising handles. When exercising the legs, the therapy grips are replaced by safety handle adapters.

**Optional fore arm rests, wristbands for the fore arm rests, tetra grips and wrist bands:** for more effective restraint of the lower arms, especially in cases of partial paralysis, there is a choice of various optional fore arm rests, wristbands for the fore arm rests, tetra grips and wristbands to suit the particular disability.

**Radius adjustment stepless option:** With the radius adjustment stepless option, the crank length and, therefore, extent of leg movement can be individually adapted to the patient without the use of tools.

**Calf rest option:** The pivoting and vertically adjustable calf rests help to secure and guide the lower legs, particularly in cases of partial paralysis.

**Foot fixing - self operated option:** The foot fixings are useful for patients to be able to quickly and easily fix their own feet into the foot rests.

**Optional tip-up protection, self-operated-tip-up protection, self-operated:** the tip-up protection, self-operated is simply hooked onto the wheelchair and tightens itself automatically. It prevents the wheelchair from tipping over when the patient is exercising with the THERA-vital.

**Optional pulse monitor with breast belt:** the Cardio-Pulse set (pulse monitor with breast belt) detects the patient's precise heart rate with the aid of a breast belt. The pulse is displayed on the large colour display.

**THERA-assist option:** the THERA-assist software is an exercising program for recording and documenting exercising data. A USB memory stick is used to transfer the data from the exerciser to the computer on which the software is installed.

## Intended Use

The THERA-vital is suitable for use in the home, in hospitals, care institutions and medical practices. It is intended to improve the mobility of persons suffering impaired mobility as a result of accidents, surgical operations, or general conditions affecting the mobility of the upper and lower limbs.

The THERA-vital can be used both as a leg exerciser and as an upper torso trainer. It is suitable for active-assisted exercise (where the patient's limbs are moved by the machine) as well as active therapeutic exercise (where the patient provides the motive power). The transition from active-assisted to active exercising is gradual. Operating both upper body and leg exerciser at the same time is not possible.

It is not possible to give specific details of the use of the THERA-vital in the case of the various possible medical conditions or to provide specific exercise programmes. The possible settings depend on the individual condition (age, build, physical constitution, fitness, ...) of the patient. Consultation with a doctor is advisable.

**The THERA-vital is a therapeutic exerciser and not a medical instrument for diagnostic purposes.**

## Exerciser Components and Controls

- 1 Control pad with display
- 2 Safety handle
- 3 Star handle for adjusting height of safety handle/upper torso trainer **H**
- 4 Footrest
- 5 Crank with two length settings
- 6 Exerciser stand, extending and height compensating **V**
- 7 Transport wheels **A**
- 8 Inclination adjuster for handlebar/upper body exerciser **I**
- 9 Transmission casing
- 10 Safety handle for upper section
- 11 Handlebar support post
- 12 Control pad inclination adjuster
- 13 Upper torso trainer
- 14 Therapy grips **K**
- 15 Crank on upper torso trainer
- 16 Star handle for adjusting reach of upper torso trainer **J**
- 17 Calf rests **F**
- 18 Crank with radius adjustment stepless **C**)
- 19 Tip-up protection **D**
- 20 Tip-up protection, self operated **E**
- 21 Foot fixing **G**
- 22 Safety handle adapter **L**
- 23 Wristband **M**
- 24 Fore arm rest for therapy grips **N**
- 25 Fore arm rests with grip bolts **O**
- 26 Fore arm rests with bar shaped grips **P** and **Q**
- 27 Wristband for the fore arm rests **R**
- 28 Tetra grips **S**
- 29 Arrester for crank **T**
- 30 Ear clip for heart rate measurement **U**
- 31 Cardio-Pulse set **U**

## Interfaces

- 32 USB interface for control pad **X**:  
A memory stick is plugged in here to update the software or to copy data for transfer to the THERA-assist program.
- 33 Fused power connection **B**:  
For connecting the power cord.
- 34 PC interface (for servicing only) **B**:  
This interface is required for commissioning and servicing.
- 35 Interface for Ear clip/Cardio-Pulse set **W**:  
This is where the optical heart rate sensor with earlobe clip **30** or Cardio-Pulse set **31** is connected.

Use only genuine medica accessories.

## Safety

- ❑ The THERA-vital is intended for supervised use only. Therefore, patients must not use the THERA-vital without being supervised by an assistant.
- ❑ Before using the exerciser for the first time or if it has been transported some distance, leave it to stand for approximately 1 hour at room temperature.
- ❑ Only move the THERA-vital on its wheels on a firm and level surface (see page 42).
- ❑ Always place the THERA-vital on a level and non-slip surface so as to ensure maximum stability. Make sure that the apparatus foot (6) is properly tightened (page 44).
- ❑ The THERA-vital may only be connected to a power supply matching the specifications on the rating plate. Only connect the exerciser to a properly earthed power socket.
- ❑ The power cord must be routed in such a way that it does not represent a tripping hazard, that it cannot become caught up in the cranks and that it cannot be damaged by other equipment. Never use the exerciser with a damaged power cord. Use only the original power cord supplied with the exerciser.
- ❑ To prevent the RISK of electric shock, only connect the exerciser to an earthed power supply.
- ❑ To prevent electric shock, the THERA-vital must never be used in wet, damp or very hot conditions.
- ❑ Have your doctor/therapist or supplier show you how to operate the exerciser before using it for the first time.
- ❑ Do not wear shoes that have laces when exercising on the THERA-vital. The shoe laces could become caught in the cranks and cause injury as a result.
- ❑ Always wear close-fitting clothing when using the exerciser.
- ❑ If the THERA-vital is to be used by a patient in a wheelchair, and if it is not possible to be absolutely certain that the wheelchair will not tip over backwards (e.g. with spasm or extremely active patients), a wheelchair spasm must be used.
- ❑ The THERA-vital is designed to be used in a sitting position only. Never stand on the foot rests 4 with the full weight of your body.
- ❑ Before commencing a leg exercising session, rotate the footrests by hand (without the motor) through one full revolution with the patient's feet in position and the legs strapped in to ensure that the cranks can move freely, that the crank length is correctly set and that movement of the legs is not restricted (legs cannot become trapped). Perform a similar check before using the upper torso trainer.
- ❑ Never reach into or attempt to take hold of moving parts of the exerciser (e.g. cranks, radius adjustment, radius adjustment, calf rests, fore arm rests, ...) when it is running.
- ❑ Before starting an upper torso trainer session, extend the apparatus foot at least 10 cm (see page 44).
- ❑ When exercising with the THERA-vital make certain that your position is physiologically correct. Seek advice from a doctor/therapist.
- ❑ In cases of partial paralysis, spasticity or instability of the legs, only use the exerciser with the optional calf rests.
- ❑ In the case of muscle contractures or limited joint movement, the radius adjustment stepless option must be used. The radius adjustment must not be altered when the motor is running. Do not adjust the radius adjustment when the patient's feet are in the foot rests. Tighten the locking screw (C) and check that it is tight before every exercising session.
- ❑ Start with an exercising session of no longer than 15 minutes. You can gradually increase the length of subsequent exercising sessions to suit the patient's individual capabilities.
- ❑ For the purposes of loosening up, always start a THERA-vital exercising session with a active-assisted exercising phase powered by the motor.
- ❑ The greater the distance between the THERA-vital and the patient's seat/wheelchair, the more the knee, hip and arm joints are extended. Therefore, you should start at a close distance. Avoid over-extending the joints.

- ❑ Before starting a leg exercising session, adjust the height of the safety handle/upper torso trainer so that the thighs can not strike any part of the exerciser. Make sure that the height adjustment (star handle, **3**) is fully tightened.
- ❑ The handlebar/upper body exerciser must not be used as a support for standing – the exerciser could tip over.
- ❑ The feet must not be placed in the footrests during an upper body exercising session.
- ❑ The adjustable-length cranks (**G**) for the upper-body exerciser **13** must not be used without additional protection against shearing.
- ❑ Make sure that the fore arm rests **24** are securely attached to the therapy grips (**N**) during an upper torso trainer session. Make sure that the arms cannot come out of the fore arm rests **24**.
- ❑ If the patient displays symptoms of illness during or after an exercising session, consult a doctor immediately.
- ❑ Children must always be constantly supervised when exercising on the THERA-vital. Never leave the THERA-vital unsupervised when children are present.
- ❑ The exerciser is only completely disconnected from the electrical power supply when the plug is removed from the power outlet socket.
- ❑ Always unplug the power cord before opening the exerciser casing. Never run the THERA-vital with the casing open.
- ❑ Repairs may only be carried out by trained specialist suppliers. If the exerciser is damaged/faulty or if it produces abnormal noises or smells, stop the exercising session immediately, unplug the power cord and contact your service engineer.
- ❑ Never use liquids containing solvents to clean the exerciser.

## Symbols Used

The symbols used in this instruction manual and, if applicable, on the exerciser are intended to draw your attention to possible dangers when using the exerciser.

You should make sure you understand the meaning of the symbols so that you can act accordingly and thereby use the exerciser more effectively and safely.

**Danger warning:** read the Instructions for Use and follow the directions given.

**Danger of moving parts:** do not reach into the cranks during an exercising session.

**Alternating current:** the THERA-vital uses an alternating current power supply.

**Earth lead connection**

## Notes on Exercise Therapy

### General observations

The object of therapy depends on the medical condition of the patient and may be aimed at maintaining mobility (preventative therapy) or at rehabilitation after an operation or injury, i.e. at regaining mobility and strength.

We recommend that the THERA-live is used as part of a planned exercise programme prescribed by a doctor or therapist.

Begin exercise sessions slowly and then increase the level of intensity gradually according to the user's physical capabilities, being particularly careful to avoid over-exertion.

### Uses

The THERA-vital can be used as a leg exerciser and as an upper torso trainer:

- Using the THERA-vital as a leg exerciser:** The patient sits on a chair or wheel chair. The THERA-vital stands on the floor in front of the patient. If possible the wheelchair footrests should be removed. The feet are held in the foot rests **4** or the calf rests **17** (optional feature), e.g. by means of the foot fixings **21** (optional feature). The exerciser is then started from the control pad **1**.
- Using the THERA-vital as an upper torso trainer:** The patient sits on a chair or wheel chair. The THERA-vital stands on the floor in front of the patient. The patient's feet are placed on the floor or the wheelchair's footrests. The upper torso trainer is adjusted to the correct height and reach.

The arms are strapped to the therapy grips **14**, wrist bands **23**, fore arm rests **24** to **27** (optional feature) or tetra grips **28** (optional feature). The exerciser is then started from the control pad **1** by an assistant.

The THERA-vital is specifically recommended for the following conditions:

- Impaired mobility and ability to walk
- General lack of mobility
- Multiple sclerosis
- Rheumatism
- Cardio-vascular problems
- Arthritis and osteoarthritis
- Parkinson's disease
- Paraplegia and tetraplegia
- Strokes and hemiplegia
- Dialysis patients
- Muscle wastage
- Persistent vegetative state/appallic syndrome

### Correct posture

Always make sure that the patient's posture when exercising with the THERA-vital is in keeping with the purpose of the therapy.

The greater the distance between the THERA-vital and the patient's chair/wheelchair, the more the knee and hip or elbow and shoulder joints are stretched. Therefore, you should start the exercise programme with the exerciser positioned close to the patient's

chair/wheelchair. In that way you will avoid overstretching the joints or damaging muscles/tendons/ligaments.

It is also important to ensure that the patient's sitting position in the chair/wheelchair is as upright as possible.

### Exercise programme planning

The frequency and duration of exercise sessions on the THERA-vital should normally be individually planned and prescribed by a doctor or therapist. Therefore, only general guidance on exercise programme planning can be given at this point.

Regular exercising with the THERA-live is extremely important if improvements in mobility, and particularly in strength and endurance, are to be achieved. Short but frequent sessions are better than long strenuous ones.

You should therefore start with sessions of no more than 15 minutes of continuous exercise. Always start with a period of gentle active-assisted exercise and then progress to light, active exercising with a low resistance setting. The length of the session, the speed, the amount of active exercising and the resistance can be gradually increased a small amount at a time.

You can also schedule several exercise sessions a day – always providing that no negative symptoms of illness occur and that the physical capabilities of the patient are not exceeded. The intensity of the therapy is correct if strength, endurance and mobility gradually improve and the patient feels well.

## Active and active-assisted exercising

**Active-assisted** exercising means that the patient's limbs are moved by the exerciser without any effort on the part of the patient. The motive power is provided by the exerciser.

When first switched on, the THERA-vital always starts with a active-assisted exercising phase with a preset resistance and speed. When the patient starts actively exercising with the same motion, the THERA-vital automatically changes over to active exercising mode.

**Active** exercising means that the patient uses his/her own efforts and energy to turn the machine's cranks against the set resistance level (see the section "Resistance/Power presetting" on page 48).

There are several different ways in which the patient can set up the machine for active exercising:

**Set resistance:** the patient exercises against a set resistance level (see page 48).

**Set power:** the patient works at a constant power output level (measured in watts). The resistance is automatically adjusted according to the pedalling speed (see page 48).

**Heart rate limit:** the patient exercises up to a level of intensity defined by a heart rate limit. If the patient's heart rate rises above the set limit, the resistance is automatically lowered.

## Muscle hypertonia/spasms

Muscle hypertonia can be effectively reduced by gentle and even movement on the THERA-vital.

The important consideration is that the distance between the THERA-vital and the patient's chair/wheelchair is kept as small as possible to begin with. If this is combined with an upright sitting position and a short crank length (see "Radius adjustment" on page 43), a reduction in the muscle hypertonia can generally be observed within a short time. Experience and scientific investigations have shown that this effect can last for several hours.

In consultation with a doctor or therapist, an exercise programme should be drawn up so as to achieve the correct "degree" of hypertonia reduction. If a spasm occurs while exercising, the electronic circuitry of the THERA-vital detects it and immediately stops the motor – assuming the anti-spasm function has been activated (factory setting).

The type of spasm suffered by most patients is an "extension spasm". On occurrence of a extension spasm, the THERA-vital stops and then continues to rotate **forwards** after a short delay. That forward rotation stretches the muscles and allows the spasm to recede (for details of settings, see page 56).

In the event of the less common type of spasm, the "flexion spasm", the THERA-vital rotates **backwards** after a spasm occurs (for details of settings see page 56).

With the **Automatic Reversal** setting (see page 56), the THERA-vital reverses the direction of rotation after a spasm has occurred.

## Transporting and Commissioning

### Before using for the first time

Remove your THERA-vital from its packaging and check for any damage suffered in transit. Also check that the power supply voltage specified on the rating plate matches that of your mains power supply.

If discover any damage or if the power supply voltage is incorrect, you should contact your supplier immediately.

### Transporting

**A** The THERA-vital is fitted as standard with transport wheels **7**. To move the exerciser, tip it as shown and push it in front of you or pull it behind you.

You can also pull the THERA-vital closer to you when sitting in a chair/wheelchair so as to obtain the best distance from the machine when exercising. Tip the THERA-vital towards your chair/wheelchair so that it is balancing on its wheels **7**. Then pull it towards you until it is at the right distance.

The wheels are not suitable for moving the exerciser over uneven surfaces or unsurfaced ground.

The exerciser must be carried up or down stairs or steps. This requires at least two persons.

### **B** Connecting to the power supply

Connect the THERA-vital to the power outlet socket using the power cord. The THERA-vital is then ready for use. You can start it from the control pad.

### Switching the exerciser on/off

If the screen is **not** illuminated, the THERA-vital is in stand-by mode. Press any button to switch the THERA-vital to normal operating mode. The screen will light up and, following a self-test lasting about 5 sec., the THERA-vital is ready for use. By pressing the green “START” button on the control pad, you can start the THERA-vital and begin exercising.

Pressing the “STOP” button on the control pad **1** during an exercising session pauses the exerciser. Pressing the button a second time ends the exercising session and the session analysis is displayed.

Pressing and holding the “STOP” button (for at least 5 sec.) sets the exerciser to stand-by mode – the screen illumination switches off

## Operation

### Foot rests

The foot rests **4** are designed for exercising with or without shoes.

The feet can be securely fixed in the foot-rests **4** with the Velcro straps provided. To do so, pull the Velcro strap over the top of the foot to the outer side of the foot rests **4** and press it against the Velcro strip.

Before every exercising session, check that your feet are correctly positioned in the foot-rests **4** and are securely held by the Velcro straps.

If you are wearing shoes, make sure that they are not open at the back because otherwise your feet could slip out of the back of the shoes and over the top of the heel retainers.

Do not wear shoes that have laces when exercising on the THERA-vital. The shoe laces could become caught in the cranks and cause injury as a result.

### Radius adjustment 2-step

The THERA-vital is fitted as standard with a radius adjustment 2-step. Using an open-ended spanner (size 13 mm) the crank length can be set to the best position. Please note that the right-hand crank has a right-hand thread and the left-hand crank a left-hand

thread. Make sure you tighten the screws properly again after changing the crank length.

### **C** Radius adjustment stepless (optional feature)

If the exerciser has been fitted with the infinitely radius adjustment, the crank length can be infinitely adjusted. Thus the degree of movement can be adjusted individually to suit each patient's mobility.

- Switch off the THERA-vital by pressing the "STOP" button.
- Remove the feet from the foot rests.
- Undo the locking screw.
- Slide the foot rest along the crank to the desired position and **firmly tighten the locking screw**.
- Adjust the opposite crank in the same way. Make sure that both cranks are set to the same length (the cranks are marked with a scale). In certain cases, it may be beneficial for therapeutic reasons to have different crank lengths on each side (e.g. if the patient's legs are different lengths or have different degrees of mobility).

For leg-amputee patients it may be necessary to completely remove one of the foot rests. To do so, undo the locking screw far enough to be able to slide the footrest together with the sliding adjuster off the end of the crank.

### **D** Tip-up protection (optional feature)

The tip-up protection prevent the wheelchair tipping over backwards.

Push the two tip-up protections firmly under wheelchair's tilting bars.

The tip-up protections can only be used with wheelchairs that are fitted with tilting bars.

### **E** Tip-up protection, self operated (optional feature)

The tip-up protection, self operated prevents the wheelchair from tipping over when the patient is exercising with the THERA-vital.

Position the wheelchair in front of the THERA-vital as required for exercising.

Press down the red lever on the first inertia belt reel **20** and at the same time pull the belt out with the hook. Attach the hook to a permanently fixed part of the wheelchair frame. The belt will then tighten itself automatically.

Do not attach the belt to removable parts of the wheelchair such as footrests, side panels, castors, etc.

The belt must be attached to the wheelchair at an angle of at least 30°.

### **F** Calf rests (optional feature)

The calf rests help to secure the legs in cases of paralysis. They are connected to the foot rests by a pivoting joint so as not to impair movement of the ankles.

- Adjust the height of the calf rests so that the cups fit around the calves. Secure the lower legs with the Velcro straps.

**Patients who are secured to the THERA-vital by the Velcro straps must not be left unattended.**

### **G** Foot fixings (optional feature)

The fixings allow the feet to be secured in the footrests **quickly, safely and in many cases by the patient**:

- Fold the foot fixing forwards.
- Place the foot in the foot rest.
- Fold the foot fixing back again.

## **H I** Adjusting the handle bar

If you wish to exercise actively but suffer from muscle hypertonia or a lack of abdominal stability, the safety handle provides extra stability and balance. You can adjust the safety handle position to suit your individual size and sitting position.

- Undo the knurled-knob locking screw **3** on the support tube.
- Adjust the safety handle to the desired height and **firmly retighten the knurled knob 3!**

**Make sure that the handlebar support post 11 is inserted to a depth of at least 10 cm in the support tube 10 – the post has a line marking the minimum insertion depth.**

- If you wish to adjust the reach of the handlebar, you must undo the eight screws **8** a small amount using the tool supplied.
- Move the safety handle to the desired position and then retighten the screws **8**.

The safety handle position should be set so as to obtain as upright an exercising position as possible.

## **V** Adjusting the apparatus foot

The THERA-vital can be given greater stability by extending the stand **6**. It is imperative to do so when using the exerciser for upper torso training.

- Tip the THERA-vital on its rear bottom edge.
- Undo the two fixing screws using the tool supplied.
- Extend the apparatus foot taking care not to exceed the maximum extension marked on the stand (12 cm).
- Retighten the fixing screws.
- Tip the THERA-vital forwards again into its normal position.

To level the THERA-vital on an uneven floor, the rubber feet on the apparatus foot can be adjusted.

- Undo the lock-nut with an open-ended spanner (size 13 mm) and unscrew the rubber foot until the THERA-vital stands firmly.
- Retighten the lock-nut.

## **Upper torso trainer (optional feature)**

### **H J** Adjusting the upper torso trainer

The upper torso trainer for the THERA-vital can be adjusted to suit the size and sitting position of the patient.

-  The upper torso trainer should be adjusted so as to obtain as upright an exercising position as possible.

## Adjusting the height

- Undo the knurled-knob locking screw **3** on the support tube.
- Adjust the upper torso trainer to the desired height and **firmly retighten the knurled knob 3.**

**Make sure that the support post 11 of the upper torso trainer is inserted to a depth of at least 10 cm in the support tube 10 – the post has a line marking the minimum insertion depth.**

**The upper section should be adjusted to a height at which the knees cannot hit the upper torso trainer at any point.**

## Adjusting the reach

The ideal reach setting for the upper torso trainer is one at which the arms are not quite fully straightened when the exerciser handles are at their furthest position from the body.

- Undo the knurled-knob locking screw **16** on the underneath of the upper torso trainer a small amount.
- Adjust the reach of the upper torso trainer as required.
- Retighten the knurled-knob locking screw 16.**

**When using the upper torso trainer, make sure that the apparatus foot 6 is extended at least 10 cm.**

**K** **Therapy grips for upper torso trainer**

The therapy grips **14** offer a variety of holding positions for the hands. The therapy grips should only be used for upper torso training if the patient is capable of gripping them sufficiently well and securely.

The therapy grips **14** are easy to remove from the spindle on the crank **15**. To do so, press in the locating pin with your thumb and pull the therapy grips **14** off the spindle. Repeat the operation with the other second therapy grips **14**.

To refit the therapy grips **14**, press in the locating pin on the therapy grips **14** and carefully slide the therapy grip onto the spindle on the crank **15** as far as the stop. Now release the locating pin. The therapy grips **14** should perceptibly snap into position.

**L** **Safety handle adapters**

The safety handle adapters convert the upper torso trainer into an upper section with safety handle in a few simple operations.

**Fitting the safety handle adapters**

- Press in the locating pin on the therapy grips **14** and pull the therapy grip off the spindle **K**.
- Press in the locating pin on the adapter handle and slide it onto the spindle on the crank **L**.
- Release the locating pin and push the safety handle adapters as far as the stop. The locating pin should audibly snap into position

**The safety handle adapters are not designed to support the full weight of the body or for carrying the exerciser.**

**T** **Locking the cranks**

- Turn the arrester knob **29** a quarter turn and then let go of it.
- Turn the cranks to a horizontal position so that the arrester knob **29** audibly engages. It should now not be possible to turn the cranks **15**.

To unlock the cranks, pull the arrester knob downwards, turn it a quarter turn and let it go.

**J** **Adjusting the position of the control pad**

The position of the control pad on the THERA-vital upper torso trainer **13** can be adjusted to suit the size and sitting position of the patient.

- Undo the screws **12** on both sides of the upper torso trainer a small amount.
- Adjust the control pad to the desired position and firmly retighten the screws **12**.

For transportation, the control pad can be folded back/down completely.

- Completely unscrew the screws **12** on both sides of the upper torso trainer.
- Fold the control pad fully back/down.
- Refit the screws **12** to fix the control pad in position for transportation.

## **N** to **S** Fore arm rests (optional feature)

The fore arm rests are a useful option for helping patients with symptoms of paralysis, hypertonia, muscle spasms or lack of hand/ arm mobility to obtain a secure hold.

Apart from fore arm rest **N** all fore arm rests can and should be individually adjusted. Combinations of fore arm rest are also possible, as is the use of only one armrest. Adjustment by a doctor or therapist is essential. Elastic Velcro straps allow the forearms to be secured in the fore arm rests.

You can either:

- place the arm in the fore arm rests after the fore arm rests has been attached to the upper torso trainer

or

- fix the fore arm rest to the arm first and then slide the fore arm rest onto the spindle on the crank.

## **M** Wristband (optional feature)

If the patient's hand is paralysed, the wristband allows it to be quickly and easily secured to the safety handle **2**, exerciser handle **14** or special Tetra-handle **28**.

## Heart rate sensor

### **U** Opto heart rate sensor with ear clip (Ear clip for pulse measurement)

The THERA-vital is supplied as standard with an ear clip for pulse measurement. The ear clip is used to monitor the heart rate (see page 56, „Pulse limit“). Due to insufficient accuracy, the ear clip must not be used to monitor the pulse limit for patients with heart and circulation conditions. In such cases the Cardio-Pulse set (optional feature) is required.

- Plug ear clip **30** **U** into the bottom of the control pad (Interface **35**, **W**).
- Then attach the ear clip to the ear lobe after first rubbing the ear lobe to ensure there is good circulation.
- Attach the coiled cord to clothing using the cable clip.

### Cardio-Pulse set (optional feature)

The chest-band heart rate sensor (Cardio-Pulse set) **31** measures the heart rate more accurately and for that reason can be used to monitor the pulse limit of patients with cardiovascular conditions (see page 56, „Pulse limit“).

- To connect it to the THERA-vital, plug the receiver connector into the socket on the bottom of the control pad (Interface **35**, **W**) so that the flat side is facing towards the patient.

- Slightly moisten the inside of the chest-band and strap it around the chest (close to the heart) so that the transmitter is facing towards the THERA-vital.

**To prevent interference when several THERA-vital exercisers are being used with Cardio-Pulse sets, there should be a gap of at least 3 m between exercisers.**

Correct functioning of the heart rate monitor can be checked on the screen: if a heart rate figure is indicated and the heart symbol is at least partially coloured red, the signals from the heart rate monitor are sufficiently strong and regular. If the heart symbol is completely coloured red, the heart rate signals are of the optimum strength.

With an earlobe heart rate sensor, the signal quality can fluctuate considerably. The signal is dependent on the blood circulation to the ear lobe and how well the sensor is attached. Wobbling of the sensor will cause interference signals and momentary breaks in signal detection that can prevent display of a heart rate even though the heart symbol is showing solid red. If no heart rate signal is present, the heart symbol is grey.

## Control pad with colour display

The colour display always shows all important settings and exercising session data. The basic start screen view shows the following information:

**Resistance level** (see page 48): shows the exerciser resistance level on a scale from 1 - 15 or the set power level from 0 - 80 watts. You can change the setting before or during the exercising session by pressing the + and – buttons.

**Power output:** only shown when exercising actively; indicates power output in watts (Neuro mode) or energy consumption in kcal (Ortho mode).

**Pulse rate,** shows the heart rate and a heart symbol that is at least partially coloured red when the heart rate monitor is connected and supplying a good heart rate signal. When the heart symbol is coloured completely red, the heart rate signal is the optimum strength. If the heart symbol is grey, no heart rate monitor is connected.

**Anti-spasm function** (see page 48): this is active if the symbol is solid red or inactive if the symbol is white.


The button next to the symbol changes the setting.

**Leg/Upper torso trainer selector** (see page 48). The active option is shown.

Change the setting as required before starting the training session press the button next to the symbol. Only possible if an upper torso trainer is fitted.

**Stops** the exercising session (see page 49)!

Pressing once pauses the exercising session. Pressing a second time ends the exercising session and displays the analysis for the last session.

**Speed** (see page 48): when exercising actively shows the actual speed of rotation of the cranks; when exercising actively-assisted shows the preset speed, which is the speed of the motor.

Press the + and – buttons to alter the setting.

**Direction of rotation** (see page 8)

Change the setting before or during an exercising session by pressing the button next to the symbol.

**Get-in-support** (see page 49).

Activate/deactivate get-in-support mode by pressing the button next to the symbol.

If get-in-support mode is active the symbol is green.

**Starts** the exercising session with an initial speed of 10 rpm and starts recording a new session analysis (see page 49).

Status bar

Remaining session time

Distance completed

### Resistance/Power presetting

You use this control to set the difficulty or power level for the exercising session. The following settings are possible:

**Resistance (Neuro mode):** In Neuro exercising mode (see page 52) you can set the resistance level.

The current resistance level is shown on a scale from 1 to 15. The maximum resistance depends on the preset motor power (see page 55) and has a maximum level of 12 Nm. The resistance can be increased by pressing the + button. To increase the resistance one increment at a time, press and release the button.

The resistance can be reduced by pressing the – button and has a minimum setting of 1.

If you press and hold the + or – button, the setting increases/decreases continuously up/down to the maximum/minimum.

The set resistance remains constant at all speeds. A higher cadence means a higher power output.

**Power presetting in watts (Ortho/Cardio mode):** in Ortho and Cardio exercising modes (see page 52) exercising is based on a power output setting (up to a maximum of 80 watts). The power level is set by pressing the + or – button. The THERA-vital automatically modulates the resistance level according to the pedalling speed so that the set power must always be delivered.

### Speed

The speed can be increased by pressing the + button. To increase the speed one increment at a time, press and release the button. If you

press and hold the + or – button, the speed increases/decreases continuously up/down to the maximum/zero. An exercising session always starts with a speed of 10 rpm.

### Spasm

The spasm detection function on the THERA-vital protects spastic patients or

those with sensitive bones or joints from undue stresses.

Pressing the button switches spasm detection on or off. When the function is active, the symbol is shown in solid red.

If the THERA-vital detects a spasm, it stops immediately. The display indicated “SPASM”. After a pause of about 5 sec. the machine continues at a active-assisted exercising speed 5 rpm lower (see also “Spasticity program” on page “Spasm program” on page 56).

Read the notes on page 41 under the heading “Muscle hypertonia/spasms”.

### Direction of rotation

The current direction of rotation is indicated by the solid arrow. Pressing the

button reverses the direction of rotation. A change of direction is always performed with a smooth transition. All other parameters (speed, motor power, arm/leg exerciser, spasticity detection) remain the same.

### Leg/Upper torso trainer

You can select either the leg exerciser or the upper torso trainer. The symbol indicates which exerciser is currently active.

You can only switch between the leg exerciser and the upper torso trainer if you have first ended the current exercising session by pressing STOP.

While an exercising session is in progress or if there is no upper torso trainer fitted, the button is deactivated.

Before starting an upper torso trainer session, the exercising handles **14** or the armrests **24 to 28** must be fitted (see page 45) and the cranks unlocked (see page 45).


## Get-in-support

You can use the get-in-support to move each pedal to an easy position for getting

the feet in.

- Press and release the button and wait until the foot rest has moved to the easy entry position (crank at lowest position).
- Place the first foot in the foot rest that is at the bottom.
- Press the button again and wait for the other foot rest to move to the easy entry position.
- Now secure the first foot in the foot rest at the top and then place the second foot in the foot rest at the bottom.
- Press the button again and wait for the foot rests to reach the first position again. Now secure the second foot.

## Changing the session time

When an exercising session is in progress, two buttons for changing the session time are shown at the bottom left and right of the screen.

Pressing the + button increases the session time. To increase the session time in one second increments, press and release the button. Pressing and holding the + or – button increases/decreases the session time continuously up/down to the maximum/minimum.

## Start/Stop

When in **stand-by** mode, the THERA-vital is activated by pressing any button. The screen then lights up. Following a self-test lasting about 5 sec. the THERA-vital is ready for use. Pressing the START button then starts an exercising session with the basic settings.

The exerciser starts at a speed of 10 rpm.

By actively turning the cranks at a speed greater than that of the motor, you can change over to active exercising at any time.

You can **pause an exercising session in progress** by pressing the STOP button once and resume it by pressing the START button.

Pressing the STOP button twice during an ongoing exercising session ends the session. The display then shows the analysis for that exercising session (see page 57).

## THERA-assist exercising software

The THERA-vital (software version 01.04.00 or later) supports the THERA-assist software for recording and documenting exercising data as standard.

THERA-assist (Art. no. 06076-000) is available as an optional feature. A detailed description of its functions is provided in the manual supplied with the software package.

## Submenu (for experienced users)

The THERA-vital has a number of exercising programs that are based on the needs of typical patient groups. For each patient group, the exercising parameters can be set on a submenu.

The complete menu structure is shown on page 51.

### Navigating the submenu

To open the submenu press the two minus buttons simultaneously.

The submenu comprising the items **Neuro**, **Ortho**, **Cardio**, **Isokinetics** and **Technical** is displayed. The function of the buttons changes as well:

**Current submenu item**

**Current selection (inverted colours)**

#### Up/More

This button scrolls the display up. In the case of menu items which require entry of an amount (e.g. motor power), this button increases the setting.

#### Down/Less

This button scrolls the display down. The menu item selected is highlighted by inverted colours. In the case of menu items which require entry of an amount (e.g. motor power), this button decreases the


#### Cancel/Discard changes

This button takes you back up to the next higher menu level.

Any changes previously saved by pressing ✓ are retained, those not saved by pressing ✗ are lost.

#### Save changes/Continue

This button saves your settings. If the menu item selected has further options, this button takes you to the next menu option.

## Menu structure

The parameters of the exercising programmes and systems settings are set out within the following menu structure:

1st Level	2nd Level	Possible Settings, Functions	Factory setting	See page
<b>Neuro</b>				
	Training time	1 to 60 minutes	15 min	52
	Forward/backward automatic	On/Off	Off	52
	Biofeedback	Symmetry bars/road, Active-assisted, Off	Symmetry bars	52
	Motor power	2 to 22 Nm	15 Nm	55
	Spasm detection	High/Medium/Low sensitivity	Medium	56
	Spasm program	Forwards/Backwards/Opposite	Forwards	56
	Automatic speed adaption	On/Off	Off	52
	Pulse limit	60 to 160	100 bpm	56
<b>Ortho</b>				
	Training time	1 to 60 minutes	15 min	52
	Biofeedback	Symmetry training bars/road, Off	Symmetry bars	52
	Motor power	2 to 22 Nm	15 Nm	55
	Muscle support	Off/On	On	56
	Pulse limit	60 to 160	100 bpm	56
<b>Cardio</b>				
	Training time	1 to 60 minutes	15 min	52
	Biofeedback	Symmetry training bars/road, Active-assisted, Off	Symmetry bars	52
	Motor power	2 to 22 Nm	15 Nm	55
	Pulse limit	60 to 160	100 bpm	56
<b>Isokinetics</b>				
	Training time	1 to 60 minutes	15 min	52
	Biofeedback	Symmetry training bars/road, Active-assisted, Off	Symmetry bars	52
	Motor power	2 to 22 Nm	15 Nm	55
	Spasm detection	High/Medium/Low sensitivity	Medium	56
	Spasm program	Forwards/Backwards/Opposite	Forwards	56
	Pulse limit	60 to 160	100 bpm	56
<b>Technical</b>				
	Microphone sensitivity	0 (off), 1 to 10	5	56
	Language			
	Contrast	1 to 40	15	
	Service settings	Security code	Four digits (0 to 9)	–

## Description of exerciser functions on submenu

The first menu level comprises the following items:

- Neuro
- Ortho
- Cardio
- Isokinetics
- Technical

### Neuro (basic setting):

Exercising programmes for patients with neurological conditions, e.g. hemiplegia, paraplegia, MS, Parkinson's disease... The method of exercising is set as **Set resistance**. The resistance level set remains constant whatever the pedalling speed. Higher pedalling speeds mean greater power output.

### Ortho:

Exercising programmes for orthopaedic patients, e.g. following a knee or hip operation. In this exercising mode, the method of exercising is preset to **Set power**. The patient is given a set power level at which to exercise. The THERA-vital modulates the resistance automatically according to the pedalling speed so that the amount of power required is always the same.

### Cardio:

Exercising programmes for cardio-vascular patients whose heart rate must be kept below a certain limit. In Cardio exercising mode, a **heart rate limit** is set. The THERA-vital reduces the resistance as soon as the patient's heart rate reaches that limit.

### Isokinetics:

This exercising program is suitable for patients with co-ordination difficulties. A **required speed of rotation** is set that is then held constant by the THERA-vital throughout the entire exercising session. The patient can thus make a physical effort by actively pushing on the pedals without having to co-ordinate the effort with the speed of rotation.

As soon as the pedalling rate drops below the required speed, the resistance is reduced. If the pedalling rate exceeds the set speed, the resistance is increased.

### Technical:

The menu option Technical allows you to make basic settings for options such as microphone sensitivity, language etc. Entry of a security code also allows access to technical service functions.

### Exercising time

The length of the exercising session can be set to between 1 and 60 minutes in all exercising modes. The default setting in each case is 15 minutes.

### Forward/backward automatic

When the Automatic direction reversal function is active, the exerciser automatically changes the direction of rotation every four minutes **when the patient is exercising actively-assisted**.

When the patient is exercising actively the direction of rotation is not reversed.

### Biofeedback leg exerciser

The purpose of the Biofeedback function is to provide the patient with feedback about the progress of the session.


The top part of the screen shows the following information:

- Current resistance/power setting
- If a heart rate sensor is connected, the current heart rate
- The patient's current power output
- The current speed of rotation of the cranks
- The distance completed
- The remaining exercising time.

### Symmetry training bars/road:


These two Biofeedback diagrams show how symmetrically **active exercising** is being performed, i.e. how great the power difference between the right and left legs for each direction of rotation (forwards and backwards, see page 48).

The bar diagram shows the distribution of power output between the two sides. The total of the two bars always adds up to 100 %.

With the road diagram, a bicycle moves to the left or right on the road according to distribution of power output between the left and right sides. In addition, the pedalling speed is indicated by the speed at which the dotted line down the middle of the road moves.

At irregular intervals, obstacles appear on the road. If an obstacle is “avoided” by appropriate “one-sided” exercising (obstacle on the left → more activity on the right) the obstacle changes to green.

The results are corrupted if...

- the cranks are not adjusted to the same length on both side (see page 42),
- the patient’s seat is not positioned centrally relative to the exerciser,
- a leg is “pulling” as well as “pushing”.


**Active-assisted mode display:**

When the patient is **exercising actively-assisted**, the colour screen shows a bar representing the motor power currently required to move the legs.

The maximum height of the bar represents 100% of the motor power setting (see page 55).

The “START” figure shows the starting level during the first minute of exercising.

The filled section of the bar and the percentage figure at the bottom (56% in this example) show what percentage of the set motor power is currently being used to move the legs.

Changing the crank length has a substantial effect on this Biofeedback display. The percentage figures from two different exercising sessions are only comparable if the crank length and the maximum motor power (sub-menu) settings are the same for both sessions.

This data is not suitable for drawing clinical conclusions. Its purpose is to show the patient how tight or loose his/her limbs are or if muscle hypertonia has been reduced.

## Biofeedback arm exerciser

The Biofeedback arm exerciser is intended to provide the patient with feedback as to how well his/her arms are working together on the cranks, i.e. how evenly they are pulling and pushing.

The purpose of the Biofeedback function is to provide the patient with feedback about the progress of the session.


The top part of the screen shows the following information:

- Current resistance/power setting,
- If a ear clip is connected, the current heart rate
- The patient's current power output
- The current speed of rotation of the cranks
- The distance completed
- The remaining exercising time.


The Biofeedback display shows two things:

- The colour indicates how well co-ordinated the two arms are during active exercising.
- The size of the coloured area indicates the patient's level of activity.


 The aim is to keep the green area as steadily as possible at the size of the thick white line.

As long as the coloured area remains green, the two arms are working harmoniously together (pulling and pushing evenly).

A change of colour to orange indicates that the arms are no longer working together harmoniously.


A sleeping smiley means that there is too little or no activity on the part of the patient.

## Motor power

The set motor power (torque) has a greater effect on the patient's legs with a smaller crank length than with a longer crank length. The motor power setting limits the motor's

maximum torque output. It can be set to between 2 and 22 Nm. The graph below provides a guide to the motor power setting based on body weight and crank length.

① to ⑥ = Infinitely adjustable crank length settings

① = Very short crank length

⑥ = Very long crank length

m = Mass (weight) of patient

M = Motor torque (power setting)

With the two-length cranks, the shorter crank length roughly equates to setting 3 and the longer crank length to setting 6.

## Example of how to use the graph:

A patient weighing 90 kg and using crank length setting 3 (for therapeutic purposes) should exercise with a motor power setting of 15 Nm.

With a short crank length (page 42) and a high motor power setting, the patient can be subjected to very strong forces. The motor power for the upper torso trainer is approx. 30% of the power of the leg exerciser.

## Spasm detection

This setting is used to specify the threshold at which the THERA-vital sensors detect a spastic reaction on the part of the patient and stops the motor to minimise the risk of overstressing the patient.

Select the setting for patients subject to spasms as follows:

**High** for patients with very sensitive joints and tendons. At this setting even very slight resistance is enough to stop the motor.

**Medium** for patients with less sensitive responses.

**Low** for patients with pronounced spasticity. The exerciser does not stop until there is a large amount of resistance.

The longer the crank length (page 42) or the lower the motor power (page 55), the more sensitively the spasm detection responds. The shorter the crank length/the greater the motor power, the less sensitively the spasm detection responds.

### Spasm program

The function spasm program is used to set the direction in which the THERA-vital rotates after detecting a spasm.

Select from the following options:

**Forwards** if the patient tends to suffer extension spasms,

**Backwards** if the patient tends to suffer flexion spasms.

**Opposite** to have the exerciser rotate in the opposite direction to the direction it was rotating when the spasm was detected

### Automatic speed adaption

If Automatic speed adaption is set to “On” the **active-assisted exercising (motor) speed** automatically adjusts to greater effort by the patient when **exercising actively**.

Example: if the patient is exercising actively at 30 rpm when the motor speed (active-assisted exercising speed) is set to 15 rpm, after about 10 seconds, the THERA-vital adjusts the active-assisted exercising speed to 25 rpm.

If the patient subsequently changes back to active-assisted exercising, the THERA-vital operates at a speed of 25 rpm. Thus the pa-

tient individually controls the speed and obtains a active-assisted exercising speed at which he/she is comfortable.

### Pulse limit

The THERA-vital reduces the resistance level as soon as the preset heart rate limit is reached. This is to ensure that the patient always exercises at a heart rate below the limit prescribed by the doctor.

 Since all Ortho and Cardio exercising programmes a high level of active effort is expected from the patient, a long crank length is advisable (see page 42).

### Muscle support

This function assists patients with only minimal residual strength or minimal strength on one side. The residual strength is supplemented by this function so that an even pedalling action is always obtained. The function operates according to the principle of “as much assistance as necessary but as little assistance as possible”. The intention is to strengthen the patient’s own efforts at all times.

In the case of leg-amputee patients, the muscle support function performs the work of the inactive side. In such cases, the muscle support function operates like an “electronic flywheel”.

### Microphone sensitivity

The THERA-vital can be stopped by a long, loud shout (safety voice deactivation). This function is particularly useful for patients who use the upper torso trainer with their arms strapped to the fore arm rests and cannot, therefore, use the control pad.

The menu option Technical/Microphone sensitivity provides the facility for adjusting the sensitivity of the microphone fitted in the control pad or switching it off altogether (indicated as  on the status bar at the bottom of the screen; the status bar only appears when the exerciser is in “Stop mode”; it is hidden while an exercising session is in progress.

Set the sensitivity as follows:

**0 (Off)** if the ambient noise is so loud that it causes the exerciser to switch off

**1 to 3** if the ambient noise is very loud

**4 to 6** if the ambient noise is of a normal level

**7 to 10** if the ambient noise is very quiet

Check that the setting is correct by carrying out a practical test.

Afterwards press the START button to continue the exercising session.

### Contrast

The contrast of the colour screen is set at the factory. If necessary, it can be adjusted by selecting the menu option Technical/Contrast.

## Keypad lock

The keypad lock can be used to prevent patients changing settings on the exerciser during an exercising session. When the keypad lock has been activated, there is a period of 30 seconds after starting an exercising session during which settings can still be changed.

- ❑ Press and hold the STOP button. Then simultaneously press the left plus button.

Successful activation of the keypad lock is indicated by the appearance of the status bar showing a lock symbol . The status bar at the bottom of the screen only appears when the exerciser is in "Stop mode"; it is hidden while an exercising session is in progress.

- ❑ To cancel the keypad lock, press the same key combination again.

**IMPORTANT: the START and STOP buttons are always active.**

## Analysis of Session Data

At the end of an exercising session, an analysis of the session data is displayed. The data is only calculated if the session has lasted for at least one minute and was started by pressing the START button.

Report	
Duration:	07:12 min
Activity:	~ 5 %
Activity left:	~ 64 %
Activity right:	~ 36 %
Distance:	0,8 km
Spasms:	3
Calories:	~ 0,0 KCal
Motor power begin:	~ 5,8 Nm
Motor power end:	~ 1,8 Nm

The following data is displayed:

- ❑ **Duration:** total length of exercising session (in minutes) excluding pauses
  - ❑ **Activity** (patient's own efforts): proportion of the total time (in %) during which the patient was exercising by his/her own efforts.
  - ❑ **Activity left:** proportion of activity for the left leg (in %).
  - ❑ **Activity right:** proportion of activity for the right leg (in %).
  - ❑ **Distance:** total distance completed (in km).
  - ❑ **Spasms:** number of spasms detected during the exercising session.
  - ❑ **Calories:** number of calories consumed (in kcal).
  - ❑ **Motor power begin:** resistance (muscle hypertonia) at the start of the exercising session (during the first 60 seconds)
  - ❑ **Motor power end:** resistance (muscle hypertonia) at the end of the exercising session
-  Certain figures (active proportion, activity left and right, calorie consumption) are only displayed if the patient has exercised actively for at least 1 minute.

Pressing the STOP or START button takes you back to the basic display.

## Activating Basic Exercising Settings

When the exerciser leaves the factory it is programmed with the basic settings as specified in the column Factory Setting on page 51.

If any settings are changed on the submenu, those changes are retained even if the power is switched off or disconnected. That means that the same settings can be used for the next exercising session.

If you wish to reset the THERA-vital to the factory settings, proceed as follows:

- Press and hold the STOP button. Then press the two minus buttons simultaneously.

The exerciser then reverts to the factory settings as specified on page 51. The settings on the Technical menu are not reset.

## Technical Data

	THERA-vital 230 V	THERA-vital 115 V	THERA-vital 100 V
Basic exerciser			
Length	74 cm		
Width	46 cm		
Height	100-110 cm		
Weight	32 kg		
With upper torso exerciser			
Length	80-100 cm		
Width	46 cm		
Height	120-130 cm		
Weight	45 kg		
Crank length			
Fixed	75 mm/ 110 mm		
Variable	65-115 mm		
Speed range	0-60 rpm		
Torque range	Approx. 2-22 Nm		
Power supply	230 V~, 50/60 Hz	115 V~, 50/60 Hz	100 V~, 50/60 Hz
Power consumption	130 VA		
Fuse	2 × 1.0 A, slow	2 × 1.6 A, slow	
Materials used	Steel, polystyrene, polyurethane, ABS and others		
Safety class	I		
Safety rating	Type B 		
Enclosure rating	I PXO		
Noise emission	Lpa ≤ 70 dB (A)		
Noise emission standard	To DIN 45635-19-01-KL2		
Ambient conditions for operation	10 °C to 35 °C 0 to 90 % Rh 970 to 1030 hPa	50 °F to 95 °F 0 to 90 % Rh 970 to 1030 hPa	
Ambient conditions for storage/transport	-30 °C to 65 °C 0 to 90 % Rh 970 to 1030 hPa	-22 °F to 149 °F 0 to 90 % Rh 970 to 1030 hPa	

## Cleaning and Disinfection

**Always make certain that you unplug the THERA-vital from the power socket before cleaning.**

Clean the surface of your THERA-vital with a soft, moist cloth.

Never use cleaning agents that are abrasive, corrosive or contain solvents, and when cleaning be careful of foils and stickers.

You can disinfect the surfaces of the THERA-vital by wiping down with a conventional disinfectant.

The machine requires no servicing.

## **B** Replacing a fuse

- Unplug the power cord from the machine
- Using a small screwdriver, remove the power supply fuse cover **33**.
- Use the screwdriver to remove both fuse holders.
- Remove blown fuses.
- Insert fuses of the correct rating as per the rating plate in the fuse holders.
- Refit fuse holder and fuses.
- Replace power supply fuse cover.

## Guarantee Terms

The guarantee period for the THERA-vital is 24 months from the date of delivery or the date of the invoice.

medica Medizintechnik GmbH undertakes to replace defective parts of the machine free of charge within that period or else to repair the machine free of charge at the factory or have it repaired by an authorised contracted agent of medica Medizintechnik GmbH.

Wearing parts are excluded from the guarantee.

Guarantee claims will not be accepted if the machine has been opened, had third-party components fitted to it or been repaired by persons not authorised by medica Medizintechnik GmbH. Damage caused by use for purposes other than that for which the machine is intended is similarly excluded from the guarantee.

## Recycling

The THERA-vital is a high-quality and durable product that has a long service, is environmentally friendly and recyclable. Most components can be recycled as used metal. Plastic components are made of polystyrene, ABS or polyurethane. The electrical and electronic components should be disposed of as electronic scrap.

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hereby declare that the product THERA-vital conforms to the following standards:

- Medicines Directive 93/42/EEC
- EMC Directive 89/336/EEC

Hochdorf, 1/ 10/ 2005

Peter Kopf

## Tips and Tricks

Problem	Suggested solution	See also
THERA-vital has too little power	Shorten crank length	Page 42
	Increase motor power one increment	Page 55
	Be certain to consult your therapist to make sure your joints and tendons are not being overstressed	
The spasm detection responds too frequently	Reduce the sensitivity of the spasm detection	Page 55
	Shorten crank length	Page 42
	Increase motor power one increment	Page 55
	Optimise sitting position/distance from exerciser	Page 40
	Be certain to consult your therapist to make sure your joints and tendons are not being overstressed	
The resistance is too high for active exercising	In Neuro exercising mode, reduce the resistance	Page 48
	In Ortho or Cardio exercising mode, reduce the power setting	Page 48
	Increase crank length	Page 42
The pulse sensor doesn't work	Make sure that the heart rate sensor is properly connected	Page 46
	Check that the ear clip is correctly attached to your ear. Is your ear warm (good circulation)? If not, briefly rub your ear lobe and/or moisten the ear clip.	Page 46
	Use clothes pegs to keep the lead in place	
	Are the transmitter (in chest band) and receiver (in the plug) facing one another?	Page 46
	Check that the chest band has been moistened and is tight enough to be in proper contact with the skin.	Page 46

<b>Problem</b>	<b>Suggested solution</b>	<b>See also</b>
The exerciser runs unevenly	Switch on muscle support	Page 56
	Set equal crank length on both sides	Page 42
	Check that sitting position is central (relative to spindle)	Page 52
	Inform service technician and have exerciser checked	
The exerciser won't start	Have you pressed START (twice if exerciser is in stand-by mode)?	Page 49
	Check power cord	Page 42
	Check power supply fuse	Page 59
	Unplug power cord and then reconnect. Then press START.	
	Select the correct exerciser (arm or leg)	Page 48
	Unlock the upper torso trainer	Page 45
Exercising session analysis is not displayed	Exercise for at least 1 minute	
Voice controlled emergency stop doesn't respond	Set higher microphone sensitivity	Page 56
Exercising settings have changed	Set new basic settings and save	Page 58
	Reset basic settings to factory settings	Page 58
Leg/Arm exerciser button has no effect	Is upper body exerciser fitted/connected?	
Symmetry training bars/street in Biofeedback display shows obviously incorrect data	Set equal crank length on both sides	Page 42
	Check sitting position is central (relative to spindle)	Page 52
	Are you "pulling" instead of "pushing" with one foot?	
	Have the THERA-vital recalibrated	
Crank doesn't stop when using get-in-support function	Press STOP button to stop accidentally started exercising session, then press and hold START button (at least 3 seconds)	
	If problem persists, contact service engineer	

## Notes on electromagnetic compatibility

### Electrical leads, lead lengths accessories

The THERA-vital may only be operated with the original power cord.

Length of the power cord: 2 m

### Warning note about accessories used

Only the original pulse sensor or Cardio-Pulse set supplied by medica may be used. Other heart rate sensors may damage the apparatus or impair its EMC characteristics.

### Warning note about siting

The apparatus or system may not be placed in the immediate proximity of or stacked together with other equipment. If operation in the immediate proximity of or stacked together with other equipment is absolutely unavoidable, the apparatus or system should be observed in order to make sure it is operating as intended in that location.

### Compliance levels

The apparatus complies with the interference immunity testing levels required by IEC 60601..

## Interference emission

Guidance and manufacturer's declaration on electromagnetic emissions		
The THERA-vital is intended for operation in the electromagnetic environment specified below. The THERA-vital owner or user should ensure that the apparatus is used in such an environment.		
Emission measurements	Compliance	Guidance on electromagnetic environment
High-frequency emissions as per CISPR 11	Group 1	The THERA-vital only uses high-frequency energy for its internal functions. Therefore, its high-frequency emissions are very low and it is unlikely that adjacent electronic equipment will suffer interference.
High-frequency emissions as per CISPR 11	Class B	The THERA-vital is intended for use in all types of facility including: - residential areas - facilities that are directly connected to a public power grid that also supplies residential buildings.
High-frequency oscillations as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker as per IEC 61000-3-3	complies	

## Interference immunity

<b>Guidance and manufacturer's declaration on electromagnetic interference immunity</b>			
The THERA-vital is intended for operation in the electromagnetic environment specified below. The THERA-vital owner or user should ensure that the apparatus is used in such an environment.			
<b>Interference immunity test</b>	<b>IEC 60601 testing level</b>	<b>Compliance level</b>	<b>Guidance on electromagnetic environment</b>
Electrostatic discharge (ESD) as per IEC 61000-4-2	± 6 kV contact discharge (indirect) ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wooden or concrete or surfaced with ceramic tiles. If the floor has a covering made of synthetic material, the relative humidity must be at least 30%.
Rapid transient electrical interference bursts as per IEC 61000-4-4	± 2 kV for power leads ± 1 kV for input and output leads	± 2 kV for power leads ± 1 kV for input and output leads	The quality of the power supply voltage should match that of a typical business or hospital environment.
Voltage surges to IEC 61000-4-5	± 1 kV out-of-phase voltage (symmetrical) ± 2 kV in-phase voltage (asymmetrical)	± 1 kV out-of-phase voltage (symmetrical) ± 2 kV in-phase voltage (asymmetrical)	The quality of the power supply voltage should match that of a typical business or hospital environment.
Voltage dips, temporary interruptions and fluctuations of the power supply voltage as per IEC 61000-4-11	< 5 % $U_T$ for 0.5 periods (> 95 % dip) 40 % $U_T$ for 5 periods (60 % dip) 70 % $U_T$ for 25 periods (30 % dip) < 5 % $U_T$ for 5 s (> 95 % dip)	< 5 % $U_T$ for 0.5 periods (> 95 % dip) 40 % $U_T$ for 5 periods (60 % dip) 70 % $U_T$ for 25 periods (30 % dip) < 5 % $U_T$ for 5 s (> 95 % dip)	The quality of the power supply voltage should match that of a typical business or hospital environment. If the THERA-vital user requires continued functioning even during power supply interruptions we recommend that the THERA-vital is supplied by an uninterruptible power supply or battery.
Magnetic field at power supply frequency (50/60 Hz) as per IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the power supply frequency should match the typical levels encountered in business and hospital environments.

NOTE:  $U_T$  is the mains alternating voltage before application of the testing level.

**Guidance and manufacturer's declaration on electromagnetic interference immunity**

The THERA-vital is intended for operation in the electromagnetic environment specified below. The THERA-vital owner or user should ensure that the apparatus is used in such an environment.

Interference immunity test	IEC 60601 testing level	Compliance level	Guidance on electromagnetic environment
Conducted high-frequency interference as per IEC 61000-4-6	3 V <sub>eff</sub> 150 kHz to 80 MHz	3 V <sub>eff</sub>	<p>Portable and mobile radio equipment should not be used in closer proximity with the THERA-vital (including power cord) than the recommended safe distance. That distance is calculated according to the equation appropriate to the transmission frequency.</p> <p><b>Recommended safe distance:</b></p> $d = [3.5/3] \sqrt{P} = 1.17 \sqrt{P}$ $d = [3.5/3] \sqrt{P} = 1.17 \sqrt{P} \text{ for 80 MHz to 800 MHz}$ $d = [7.0/3] \sqrt{P} = 2.33 \sqrt{P} \text{ for 800 MHz to 2.5 GHz}$ <p>where P is the rated output of the transmitter in watts (W) as stated by the transmitter manufacturer and d is the recommended safe distance in metres (m).</p> <p>The field strength of fixed radio transmitters at all frequencies should be below the compliance level based on on-site measurements.</p> <p>Interference is possible in the vicinity of equipment that carries the following symbol.</p> <div style="text-align: center;">  </div>
Irradiated high-frequency interference as per IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

NOTE 1: At 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2: This guidance may not be applicable in all cases. The propagation of electromagnetic energy is affected by absorption and reflection by buildings, objects and people.

- ❑ The theoretical field strength of fixed transmitters such as radio telephone and mobile agricultural radio equipment base stations, amateur radio transmitters, AM and FM radio and TV transmitters cannot be precisely determined in advance. In order to determine the electromagnetic environment in respect of fixed transmitters, a study of the location should be considered. If the measured field strength at the location of the THERA-vital exceeds the compliance levels, the THERA-vital should be observed to verify that it is functioning as intended. If unusual performance characteristics are observed, additional measures may be required such as modifying or changing the location of the THERA-vital.
- ❑ Above the frequency range of 150 kHz to 80 MHz the field strength should be less than 3 V/m.

### Recommended safe distances between portable and mobile equipment, high-frequency telecommunications equipment and the THERA-vital

The THERA-vital is intended for operation in the electromagnetic environment specified below. The THERA-vital owner or user can help to prevent electromagnetic interference by observing the minimum distances between portable and mobile high-frequency telecommunications equipment (transmitters) and the THERA-vital as specified below (according to the output power of the communication device).

Rated output of transmitter [W]	Safe distance based on transmission frequency [m]		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.17 \sqrt{P}$	$d = 1.17 \sqrt{P}$	$d = 2.33 \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30

For transmitters whose rated output is not specified in the above table, the distance can be determined using the appropriate equation in each case and in which P is the rated output of the transmitter in watts (W) as stated by the transmitter manufacturer.

**NOTE 1 :**

For calculating the recommended safe distance from transmitters within the frequency range of 80 MHz to 2.5 GHz an additional factor of 10/3 was applied. That reduces the likelihood that a mobile/portable communications device inadvertently brought into the patient area will cause interference.

**NOTE 2:**

This guidance may not be applicable in all cases. The propagation of electromagnetic energy is affected by absorption and reflection by buildings, objects and people.

**Index**

<b>A</b>		<b>C</b>		<b>E</b>
Active exercising .....	41	Calf rests .....	43	EC Declaration of Conformity .....
Active proportion .....	57	Calories		Exercise programme planning .....
Active-assisted exercising.....	41	Consumed (exercising session		Exercise therapy notes .....
Active-assisted exercising speed.....	48	analysis).....	57	Exerciser components .....
Activity.....	57	Cardio .....	51	Exerciser options .....
Left (exercising session analysis) ....	57	Cardio-Pulse set .....	46, 62	Exercising
Right (exercising session analysis) ..	57	Cleaning and disinfection.....	59	active.....
Adjusting crank length.....	42	Code number .....	51	active-assisted .....
Adjusting height of upper section.....	44	Commissioning .....	42	Method of exercising.....
Analysis of session data .....	57	Connecting to the power supply.....	42	Session data .....
Arm/Leg exercising .....	48	Control pad .....	47	Exercising handles .....
Arm/Upper torso trainer		Controls.....	37	Exercising level .....
locking/unlocking.....	45	Correct posture .....	40	Exercising speed.....
Automatic speed adaption .....	56			Exercising time.....
		<b>D</b>		Extension spasms.....
<b>B</b>		Direction of rotation.....	48	
Basic settings (factory settings) .....	51, 58	Forward/backward automatic.....	52	<b>F</b>
Biofeedback		Disinfection and cleaning .....	59	Factory settings.....
Arm exerciser.....	54	Distance		Flexion spasms .....
Leg exerciser .....	52	Distance completed .....	47	Foot fixing .....
Buttons		Total (exercising session analysis) ..	57	Foot fixings.....
Description of function .....	47, 50			Fore arm rests.....
Start .....	49			
Stop.....	49			Guarantee terms .....
<b>G</b>		Get-in-support .....	49	
General observations.....	40	Guarantee period .....	59	

<b>I</b>		<b>O</b>	
Intended use .....	36	Operation .....	42
Isokinetics .....	51, 52	Ortho .....	51
<b>K</b>		<b>P</b>	
Keypad lock .....	57	Power supply .....	58
<b>L</b>		Power supply fuse.....	58
Lef/Upper torso trainer		Pulse	
Switching between .....	48	Pulse limit.....	51, 52, 56
Leg exerciser .....	40	Pulse sensor .....	46
Levelling.....	44	Pulse rate	
Locking the cranks .....	45	Pulse rate display.....	47
<b>M</b>		<b>R</b>	
Microphone sensitivity for safety voice deactivation.....	56	Recycling .....	59
Microphone sensitivity for safety voice deactivation .....	51	Replacing a fuse .....	59
Motor power .....	55	Resistance .....	48
Muscle hypertonia (exercising session analysis).....	57	Resistance level.....	47
Muscle power assistance .....	56		
Muscle support.....	56		
<b>N</b>			
Neuro .....	51		
Notes on electromagnetic compatibility .	62		
Number of spasms (exercising session analysis).....	57		

<b>S</b>		Submenu.....	50
safety .....	56	Menu structure .....	51
Safety handle .....	44	Navigating .....	50
Safety handle adapters .....	45	Switching on/off.....	42
Sensitivity of microphone .....	56	Symmetry diagram.....	52
Sensitivity of spasm detection.....	55	<b>T</b>	
Servicing .....	59	Technical data.....	58
Session time .....	49	THERA-assist exercising software.....	49
Set		Tips and tricks.....	60, 61
Power (Ortho) .....	52	Total distance (exercising session	
Pulse limit (Cardio).....	52	analysis).....	57
Pulse limit (Cardio) (Cardio).....	56	Transporting.....	42
Resistance (Neuro).....	52	<b>U</b>	
Set power (Ortho) .....	52	Upper body exerciser	
Set resistance .....	52	locking/unlocking.....	45
Spasm.....	48	Upper body exercising .....	45
Sensitivity of spasm detection.....	55	Upper torso trainer .....	44, 45
Spasm:Spasm program .....	56	User-operated stabiliser .....	43
Spasticity		Uses.....	40
Anti-spasm function .....	47	<b>W</b>	
Speed		Wristband.....	46
Active-assisted exercising speed.....	48		
Automatic adaption .....	56		
Stabilizer			
User-operated .....	43		
Start .....	49		
Start/Stop button .....	49		
Stop.....	49		